



Monarch[™] Gel Extraction / DNA Cleanup Kits


- $\geq 6\mu\text{l}$ 超低量Elute體積
- 磨砂管面與標籤欄標記方便
- 縮短離心及incubation時間
- Buffer使用不同顏色辨識,防止操作失誤
- 最環保的設計,塑膠耗材減量44%及使用再生紙包裝
- Buffer與膜管均可分別購買




NEB特殊低殘留膜管設計，可防止鹽類及液體殘留



NEB之核酸回收率優於Q產品與M產品



NEB之核酸回收濃度優於Q產品與M產品



NEB之回收效果較Q牌更佳

Catalog #	Products	Kit Components
T1020S/L	Monarch DNA Gel Extraction Kit (50/250)	Gel Dissolving · DNA Wash Buffer · DNA Elution Buffer · DNA Cleanup Columns (5 µg)
T1030S/L	Monarch PCR & DNA Cleanup Kit (50/250)	DNA Cleanup Binding Buffer · DNA Wash Buffer · DNA Elution Buffer · DNA Cleanup Columns (5 µg)