



Quick Blunting™ and Quick Ligation™ Kits

- 可讓各種DNA片段快速形成blunt end，並進行Blunt end Ligation
- 室溫反應，Blunting反應僅需15分鐘，Ligation反應僅需5分鐘
- 可無需考慮切位互補性,直接以blunt end特性進行連接
- 適用於各式Cloning ,Library建構及所有傳統Ligation反應



室溫進行ligation反應5分鐘，即可穩定獲得良好cloning結果

Catalog #	Products	Kit Components
E0542 S/L	Quick Blunting and Quick Ligation Kit (20/100 Rx)	
E1201 S/L	Quick Blunting Kit (20/100 Rx)	Blunting Enzyme Mix, 10X Blunting Buffer, dNTP Mix
M2200 S/L	Quick Ligation Kit (30/150 Rx)	Quick T4 DNA Ligase, Quick T4 DNA Rx Buffer